de	Gebrauchsanweisung	4
en	User's Manual	13
fr	Mode d'emploi	22
it	Istruzioni per l'uso	31
es	Instrucciones de empleo	40
pt	Manual de operação	49
ru	Руководство по применению	58

MASTERLIGHT®

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User's Manual – MASTERLIGHT® Examination Lamps LED, 18W and 18W, Focus

Floor stand model with pivoting arm: 18W – with fixed focus and 18W, Focus – with focusable light spot

Dear customer, congratulations on your purchase of the new KaWe examination light. The innovative technology and design of this new generation of LED examination lights will help you to keep your professional edge.

Advantages of LED technology: a lamp life of at least 50,000 hours and hardly any noticeable heat build-up near the doctor's head or on the area of the patient's body that is being examined. The advantages of conventional illumination technology with halogen and gas-discharge light bulbs have been maintained: natural colour reproduction, exact illumination of the treatment area and easy positioning of the light.

► TABLE OF CONTENTS

1. Safety instructions	page 13	5.1 Procedures for first use and maintenance	page 16
2. Short description of the KaWe examination light	page 14	6. Technical data	page 17
3. Operating the KaWe examination light	page 14	6.1 Photometric data	page 17
3.1 Switching the light ON/OFF	page 14	6.2 Electrical data	page 17
3.2 Brightness control	page 14	6.3 Electrical installation information	page 17
3.3 Focussing (only 18W, Focus)	page 14	6.4 Weights	page 17
3.4 Positioning	page 14	6.5 Environmental conditions	page 17
4. Cleaning	page 15	6.6 Important information	page 18
4.1 Sterilisable handle	page 15	7. CE-label	page 18
4.2 Light body, protective screen and support system	page 16	8. Disposal	page 18
5. First use and maintenance	page 16	9. Electromagnetic compatibility tables	page 18

► 1. SAFETY INSTRUCTIONS

When handling the light, the instructions given in the User's Manual must be followed.

CAUTION : This device is not intended for use in areas in which explosions may occur. This lamp is classified as a class 1 device according to the medical products law.

Prior to assembly, store the KaWe examination light in its packaging in the room in which it will be used for at least 24 hours in order to compensate for transport-related temperature fluctuations.

Please read these instructions thoroughly and carefully and familiarise yourself fully with this product before using it in order to benefit from all the advantages of the lamp system and to avoid the possibility of damaging it.

Repair and maintenance work on the KaWe examination light is only to be carried out by us or by a location expressly authorised by us to perform such work.

The safety of the KaWe examination light is only guaranteed if repairs and changes are performed by the manufacturer itself or by a location that guarantees that the safety rules will be followed.

⚠ The KaWe examination light may not be modified in any way!

The manufacturer is not responsible for damage to persons or other objects if the KaWe examination light is used incorrectly or for a non-intended purpose. Disassembly of the light body from the spring arm is to be carried out in reverse order of assembly and may only be carried out after the arm has been raised to its maximum height and is in a horizontal position.

Prior to each use, check the KaWe examination light to ensure that it is in perfect working order.

Warning, external adapter! The KaWe examination light only works with the 60VA external adapter. The external adapter that is used for the examination lamp must be tested in accordance with IEC 60601-1.

ON/OFF switch: An on-site ON/OFF switch is to be installed to interrupt the input current to the device. The switch must fulfil IEC 61058-1 requirements for rated voltage peaks of 4kV.

CAUTION: During assembly of the KaWe examination light, the entire system (including ceiling mount) must be disconnected from the power supply! Later removal of the KaWe examination light from the spring arm or disassembly of the sliding contacts in the arms is **ONLY ALLOWED AFTER COMPLETELY DISCONNECTING THE LIGHT FROM THE POWER SUPPLY NETWORK.**

The electronics may become damaged if the light is not disconnected from the power supply network!

Symbol key:

	Caution!
	Important assembly instructions, useful information and user tips
	Temperature limit
	Separate disposal of electric and electronic devices
	Complies with relevant EU guidelines
	Heed the User's Manual
	Alternating current
	Bulb
	Protective earth connection
	Protection class II
	Fuse
	Manufacturer
	Series number
	Product reference number
	Non-ionizing radiation

2. SHORT DESCRIPTION OF THE KaWe EXAMINATION LIGHT

Information about intended use: The KaWe examination lights were developed to be used by doctors in hospitals or in clinics in order to illuminate areas of patient's bodies during medical examinations.


Main function: The KaWe examination light serves to provide illumination.

General product description

- This is a KaWe examination light as described in EN 60601-2-41, which as individual luminaries, are not failsafe.
- The KaWe examination light is intended to aid in treatment and diagnosis.
- The KaWe examination light is used in areas used for medical reasons (group 0, 1 and 2 according to DIN VDE 0100-710 and HD 60364-7-710).
- The lights are mounted on a stand.
- The lights are to be serviced every 2 years.

The KaWe examination light is available in the following versions: 18W (fixed focus) – with electronic brightness control and 18W, Focus – with focusable light spot and electronic brightness control.

3. OPERATING THE KaWe EXAMINATION LIGHT

3.1 Switching the light ON/OFF

The KaWe examination light is turned OFF and ON using button **1** on the control panel.

3.2 Brightness control

The KaWe examination lights are all equipped with a brightness control function. The lamp models offer brightness modulation between 50% and 100% allowing the brightness of the light to be selected as required. By pressing button **2**, the brightness can be reduced. Button **3** increases the brightness. The selected intensity is visible on the display **4**.

3.3 Focussing (only 18W, Focus)

The KaWe examination light 18W has a light spot that can be focused, meaning it can be enlarged or made smaller as the situation requires. To focus the light spot, turn the handle **5** on the light body (see image).

3.4 Positioning


Use the handle **5/6** or both handlebars **7** to position the light body as desired.

The outer handlebars are used to position the light before a procedure. The handle is used to reposition the KaWe examination light during surgical procedures.

There are two types of handles available:

- Standard handle **5**
- Sterilisable handle **6** (additional cost)
The sterilisable handle can be removed for sterilisation.

► 4. CLEANING

4.1 Sterilisable handle

For an additional charge, the KaWe examination light can be equipped with the **sterilisable handle sleeve 8**. The removable handle sleeve can be steam sterilised and must be cleaned, disinfected and sterilised before the first and each subsequent use.

To be sterilised, the handle sleeve must be removed:

- To remove the sterilisable handle sleeve **8**, press and hold down the locking mechanism **V** and pull the handle sleeve off in a downward motion.
- Put the handle sleeve **8** back on by pushing it back on with a twisting motion until the locking mechanism **V** clicks securely.

 During an examination, the handles often become unsterile. Therefore, ensure that additional replacement handles are kept at hand.

Cleaning/Disinfection and Sterilisation

Basics

Effective cleaning and disinfection are key requirements for effective sterilisation of the handle. The responsibility for the sterility of these products includes ensuring that only sufficiently-validated equipment and product-specific processes are used for cleaning/disinfection and that the validated parameters are complied with during every cycle. In addition, the hospital / clinic hygiene regulations must be observed.

Note: The national board requirements (norms and directives) for hygiene and disinfection must be followed.

Cleaning/disinfection

Cleaning and disinfection must take place immediately after use. A mechanical sterilisation procedure should be used for cleaning / disinfection. Only procedures with proven effectiveness may be used (e.g. those listed under disinfectants and disinfection procedures that are tested and recognised by Robert-Koch-Institute / DGHM) and they must always have already been previously validated. When using other procedures (such as manual cleaning), proof of the effectiveness of the procedure must be given as part of the validation. The principle proof of the suitability of the handles for efficient cleaning / disinfection was carried out using a cyclic cleaning system (Netsch-Bellmed T-600-IUDT/AN, programme 2 for small parts; code B).

Agents / disinfectants containing the following substances may not be used as these may cause changes in the material:

- High-concentration organic and inorganic acids
- Chlorinated hydrocarbons
- 2-ethoxyethanol

When cleaning and disinfecting, the following procedures must be followed:

Zone	Procedure	Time (sec.)
1	Pre-rinse, external, cold, 10 – 15°C	45
	Washing, acidic, external 35°C	120
	Draining time	10
	Re-rinse, external, approx. 80°C	*10
	Draining time	*15
	Re-rinse, external, approx. 80°C	*15
	Draining time 45	15
2	Washing, alkaline, external, 93°C	135
	Draining time	10
	Re-rinse, external, approx. 90°C	10
	Draining time	15
	Re-rinse, external, approx. 90°C	15
	Draining time	15
3	Drying, external 100 – 120° C	200
4	Drying, external 100 – 120° C	200
	Door open / close & transport (sluice discharge)	60
	Total cycle time approx.	290 ≈ 5 Min.

* During the disinfection zone (wash zone 2), the re-rinse and draining times will depend on the type of objects being washed!

Sterilisation

Only previously cleaned and disinfected handles may be sterilised. The handles are placed in and sterilised in suitable sterilisation packaging (single-use sterilisation packaging, such as wrap / paper sterilisation bags, single or double pack) in accordance with DIN EN 868 / ISO 11607 for steam sterilisation. Only the sterilisation procedure listed below may be used. Other sterilisation procedures (such as ethylene oxide, formaldehyde and low-temperature plasma sterilisation) are not permissible.

Steam sterilisation procedure

Validated in accordance with DIN EN 554/ISO 11134

Maximum sterilisation temperature 134°C

The principle proof of the suitability of the handles for effective sterilisation was carried out using a fractional vacuum process (Euroselektomat 666 by MMM Münchner Medizin Mechanik GmbH, sterilising temperature 134°C, holding time 7 min.). If other procedures (such as manual cleaning) are used, proof of the suitability and principle effectiveness of the procedure must be given as part of the validation.

Inspection/durability

Each time before they are reused, the handles should be inspected for damage and replaced if necessary.

The sterilisable handle sleeve must be discarded and replaced with a new one after 1000 sterilisation cycles or after a maximum of 2 years. A stamped number on the inside of the handle sleeve (as shown) can be used to check the year of manufacture of the sleeve (in the picture, 12 stands for the year 2012).

4.2 Light body, protective screen and support system

The surface of the light is made of a high-quality material that can be cleaned with conventional cleaning agents.

The protective screen **9** is made of a high-quality plastic. Pay attention to the following during cleaning:

- Always clean the screen **9** with a wet cloth (never clean with a dry cloth!).
- The protective screen **9** is to be wiped with an anti-static agent using a lint-free cloth.

► **5. FIRST USE AND MAINTENANCE**

The KaWe examination light is to be serviced at least every two years. This includes an electrical and mechanical inspection.

Please heed the user's manual and assembly instructions for the support systems as these may contain different maintenance timeframes.

Caution: Make sure that the height stop (if present) is horizontal before the lamp is removed from the spring-loaded arm. Please heed the user's manual and assembly instructions for the support system.

 **Note:** Before conducting any maintenance or tests on the KaWe examination light, unplug it from the electrical outlet and ensure that it will not be plugged in again until the work has been completed.

5.1 Procedure for first use and maintenance

Inspect the KaWe examination light (especially for the following):

- Paint defects
- Cracks in plastic parts
- Deformation of the support system
- Loose parts
- Check the connection between the examination lamp and support system.
- Inspect and grease the mechanism securing the arm to the stand
- Proper operation
- Electrical safety

Note: Circuit diagrams, part lists and maintenance instructions are available upon request. The replacement of parts and repair work on the lamp during treatment are not allowed. Touching parts contained within the lamp housing while at the same time touching a patient is not allowed.

► 6. TECHNICAL DATA

6.1 Photometric data

	18W, Focus	18W
Central illumination intensity at a distance of 1 meter	45,000 Lux	35,000 Lux
Light spot diameter d_{10}	122 mm	132 mm
Light spot diameter d_{50}	62 mm	66 mm
Residual illumination intensity with one shutter	0 %	0 %
Residual illumination intensity with two shutters	60 %	61 %
Residual illumination intensity at the bottom of a normed tube	100 %	100 %
Residual illumination intensity at the bottom of a normed tube and one shutter	0 %	0 %
Residual illumination intensity at the bottom of a normed tube and two shutters	60 %	61 %
Illumination depth 20 %	1750 mm	1750 mm
Illumination depth 60 %	890 mm	840 mm
Colour rendering index R_a	95	95
Colour rendering index R_g	94	94
Max. irradiation intensity in the field at a distance of 1.0 m	162 W/m ²	126 W/m ²
Max. irradiation intensity in the field at a distance of 0.8 m	201 W/m ²	148 W/m ²
Focusable light spot size	14-20 cm	17 cm
Colour temperature (Kelvin)	4300 K	4300 K
Heat build-up near the head	0.5° C	0.5° C
Electronic brightness control on the light body (Standard)	50 - 100 %	50 - 100 %
Number of LEDs	12	12
Working lifespan of the LEDs	50,000 h	50,000 h
Work area	70-140 cm	70-140 cm
Diameter of the light body	29 cm	29 cm
Height adjustment limit	121 cm	121 cm

Note: The technical data are subject to certain deviations. Due to production-related reasons, the actual values can vary slightly from the above-named values. The values for R_a can deviate by approx. $\pm 5\%$. The colour temperature values can deviate by approx. $\pm 200K$.

6.2 Electrical Data

	18W, Focus / 18W
Power input	18 W
Operating voltage	24 V DC
Amperage	0.75 A

6.3 Electrical installation information

The KaWe examination light is exposed to a current peak when it is switched on and is therefore equipped with a power adapter. The adapter has a wide range input, input voltage 100 - 240V AC, 50 - 60Hz, Output voltage 24V DC.

If a changeover relay is required for an existing on-site emergency power supply for the KaWe adapter, this relay must be ordered separately from KaWe.

Warning! This light is a Protection Class I device. To avoid the risk of electric shock, this unit must be connected to a grounded electric outlet.

6.4 Weights

KaWe examination lamp	Weight
18W	1.9 kg
18W, Focus	1.9 kg

6.5 Environmental conditions

	Operation	
	Min.	Max.
Temperature	+10° C	+30° C*
Relative humidity	30 %	75 %
Air pressure	700 hPa	1060 hPa

* for higher temperatures please contact us

	Transport/storage	
	Min.	Max.
Temperature	-10° C	+50° C
Relative humidity	20 %	90 %
Air pressure	700 hPa	1060 hPa

Instructions on the packaging

Temperature range during transport and storage 	Humidity during transport and storage RH 20% - 90%	Air pressure during transport and storage P 700hPa - 1060hPa
--	---	---

6.6 Important information

If multiple KaWe examination lights are used at once, the total irradiance may exceed $1,000 \text{ W/m}^2$ due to superposition of the luminous fields. As a result, there is a risk that the illuminated field may become very hot. The superposition of the luminous fields of multiple KaWe examination lights may cause the limit values for UV radiation ($<400 \text{ nm}$) of 10 W/m^2 to be exceeded.

The report on the factory test for electrical safety can be obtained upon request. The only requirement is the serial number of the KaWe examination light for which the protocol is desired. If other KaWe examination lights or other pieces of equipment are connected together during installation, section 16 of EN 60601-1:2013 applies and if necessary, the compliance with the requirements is to be verified. Before initial operation, the installation must be inspected in accordance with EN 62353.

When installing the KaWe examination light, the polarity is very important. Should the KaWe examination light not work after installation, during troubleshooting, the polarity can be switched on the secondary side of the power supply.

7. CE-LABEL

This KaWe examination light complies with guidelines 93/42/EEG (guidelines for medical products from the council of the European communities). The applicable standard is EN 60601-2-41. The company KaWe GmbH + Co. KG is certified according to DIN EN ISO 13485:2012 + AC:2012.

8. DISPOSAL

At the end of the product's life, the components of the KaWe examination light are to be disposed of properly. Ensure that the materials are sorted carefully according to type. The electrical circuit boards are to be disposed of at an appropriate recycling centre. The lamp housing and the rest of the components of the KaWe examination light should be disposed of according to their type of material.

9. ELECTROMAGNETIC COMPATIBILITY TABLES AND INFORMATION

These KaWe examination lights are subject to special safety measures with regard to EMC requirements and must be installed in accordance with the enclosed EMC instructions. The functionality of these KaWe examination lights can be influenced by portable and mobile HF communication devices.

The use of other accessories will result in increased emissions or decreased immunity of the unit.

For the KaWe examination light to operate as intended, it is necessary that the light not be placed directly next to or stacked together with other equipment. If the lamp must be used directly next to or while stacked with other equipment, the KaWe examination light is to be closely monitored.

Table 9.1 Guidelines and Manufacturer Declaration – Electromagnetic Emissions

Guidelines and Manufacturer Declaration – Electromagnetic Emissions		
The KaWe examination lamp is designed to be used in the types of environments listed below. The customer or user of the KaWe examination lamp is responsible for ensuring that this device is used in such an environment.		
Emissions Measurement	Conformity	Electromagnetic environment – guidelines
Harmonic oscillations according to IEC 61000-3-2	Class A	The KaWe examination lamp is designed for use in all facilities including living spaces that are directly connected to a public low-voltage power supply network that also supplies power to buildings used for domestic purposes.
Voltage fluctuations / flicker emissions according to IEC 61000-3-3	Complies	
Type-CISPR 15-1 high frequency emissions	Complies	The KaWe examination lamp is not designed to be connected to any other type of device.

Table 9.2

Guidelines and Manufacturer Declaration – Electromagnetic Interference Immunity

Guidelines and Manufacturer Declaration – Electromagnetic Interference Immunity			
The KaWe examination lamp is designed to be used in the types of environments listed below. The customer or user of the KaWe examination lamp is responsible for ensuring that this device is used in such an environment.			
Interference immunity test	IEC 60601-test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) according to IEC 61000-4-2	$\pm 6\text{ kV}$ contact discharge $\pm 8\text{ kV}$ air discharge	$\pm 6\text{ kV}$ contact discharge $\pm 8\text{ kV}$ air discharge	The flooring should be wood or concrete or be covered with ceramic tile. If the floor material is a non-conductive, synthetic material, the relative humidity of the air must be at least 30%.
Fast transient electrical disturbances/ bursts according to IEC 61000-4-4	$\pm 2\text{ kV}$ for power supply lines $\pm 1\text{ kV}$ for input and output lines	$\pm 2\text{ kV}$ for power supply lines Not applicable	The quality of the mains supply voltage should be the same as that of a typical commercial or hospital environment.
Surges according to IEC 61000-4-5	$\pm 1\text{ kV}$ differential mode voltage $\pm 2\text{ kV}$ common mode voltage	$\pm 1\text{ kV}$ differential mode voltage $\pm 2\text{ kV}$ common mode voltage	The quality of the mains supply voltage should be the same as that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines according to IEC 61000-4-11	$< 5\% U_T$ ($>95\%$ dip in the U_T) for a $\frac{1}{2}$ period $40\% U_T$ (60 % dip in the U_T) for 5 periods $70\% U_T$ (30 % dip in the U_T) for 25 periods $< 5\% U_T$ ($>95\%$ dip in the U_T) for 5 seconds	$< 5\% U_T$ ($>95\%$ dip in the U_T) for a $\frac{1}{2}$ period $40\% U_T$ (60 % dip in the U_T) for 5 periods $70\% U_T$ (30 % dip in the U_T) for 25 periods $< 5\% U_T$ ($>95\%$ dip in the U_T) for 5 seconds	The quality of the mains supply voltage should be the same as that of a typical commercial or hospital environment. If the user of the KaWe examination lamp requires continued operation during power interruptions, it is recommended that the KaWe examination lamp be powered from an uninterruptable power supply or a battery.
Magnetic field with a power frequency of (50/60 Hz) according to IEC 61000-4-8	3 A/m	30 A/m	Magnetic fields with their line frequencies should be the same as those typical of commercial and hospital environments.
Note: U_T is the mains voltage prior to application of the test level.			

Table 9.3

Guidelines and Manufacturer Declaration – Electromagnetic Interference Immunity

Guidelines and Manufacturer Declaration – Electromagnetic Interference Immunity			
The KaWe examination lamp is designed to be used in the types of environments listed below. The customer or user of the KaWe examination lamp is responsible for ensuring that this device is used in such an environment.			
Interference immunitytest	IEC 60601-test level	Compliance level	Electromagnetic environment – guidelines
Conducted HF-interference according to IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V	Portable and mobile HF communications equipment should be used no closer to the KaWe examination lamp (including its power cords) than the recommended separation distance. The separation distance is calculated using various equations depending on the transmission frequency. Recommended separation distance: $d = 1.17\sqrt{P}$
Radiated HF-interference according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.17\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.34\sqrt{P}$ for 800 MHz to 2.5 GHz “P” is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. “d” is the recommended separation distance in meters (m). Field strengths from fixed transmitters, as determined by an electromagnetic survey of site “a” are less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1: For 80 MHz and 800 MHz, the larger of the two values applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted with accuracy in theory. It is recommended that in order to assess the electromagnetic environment caused by fixed HF transmitters, an electromagnetic site survey should be conducted. If the measured field strength at the location at which the KaWe examination lamp is used exceeds the applicable compliance level stated above, the equipment should be checked at each of its locations of use in order to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the KaWe examination lamp. b: For frequencies ranging from 150 kHz to 80 MHz, the field strength is less than 3 V/m.			

Table 9.4

Recommended Separation Distances between Portable and Mobile HF Communications Equipment and the device or system

Recommended Separation Distances Between Portable and Mobile HF Communications Equipment and the KaWe examination lamp			
The KaWe examination lamp is intended for use in an electromagnetic environment in which radiated HF disturbances are controlled. The customer or the user of the equipment can help prevent electromagnetic interference by ensuring that the minimum distance (shown below) between communications equipment (transmitters) and the KaWe examination lamp is maintained. The minimum distance is dependent on the maximum output power and the frequency of the communications equipment.			
	Separation Distance According to Transmitter Frequencies (m)		
Rated Power Output of Transmitter (W)	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.34\sqrt{P}$
0.01	0.12	0.12	0.23
0.10	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
For transmitters with a maximum rated output not specified in the above table, the recommended separation distance (d) can be determined in meters (m) using the equation for each column, where the maximum rated power (P) of the transmitter is in watts (W) as given by the transmitter manufacturer.			
Note 1: For 80 MHz and 800 MHz, the larger of the two values applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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KIRCHNER & WILHELM GmbH + Co. KG